



610.756.7879

[Home](#) » [Products](#) » [Merit Bergh Cilia Forceps](#)

MERIT BERGH CILIA FORCEPS

Product Number: 98-155

3-1/2" Merit Bergh Cilia Forceps Horizontal Serrations



Additional Information

Surname	Bergh
Length	3-1/2"
Tip Configuration	Blunt
Instrument Type	Tissue Forceps
Curvature	Straight
Working Surface Style	Horizontally Serrated
Handle	Thumb
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Physician Grade (Pakistan)

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7879 • Fax: 610.430.3941 • Service@SklarCorp.com