



610.756.7884


[Home](#) » [Products](#) » [Sklar Edge TC Metzenbaum-Lahey Dissecting Scissors - 6"](#)

## SKLAR EDGE TC METZENBAUM-LAHEY DISSECTING SCISSORS - 6"

Product Number: 16-1905

6" Sklar Edge Tungsten Carbide Metzenbaum-Lahey Dissecting Scissors, Curved, Delicate

Sklar Edge™ Tungsten Carbide Metzenbaum-Lahey Dissecting Scissors are most commonly used for dissecting and cutting tissue. Metzenbaum-Lahey scissors tend to be shorter than traditional Metzenbaum scissors. Sklar Edge™ Tungsten Carbide Scissors provide an exceptionally durable and hard cutting edge. These scissors are designed to stay sharp for a significantly longer period of time than standard scissors. Sklar permanently bonds tungsten carbide inserts to the blades of the scissors which maximizes control and precision. Sklar Edge™ scissors can be identified by two gold coated finger rings. This product has curved, smooth, blunt tips, and an overall length of 6 inches. It also has a delicate pattern.



### Additional Information

Surname	Metzenbaum-Lahey
Length	6"
Tip Configuration	Blunt/Blunt
Instrument Type	Dissecting Scissors
Curvature	Curved
Working Surface Style	Smooth
Pattern	Delicate
Handle	Finger Rings / Dual Gold
Material	Tungsten Carbide / Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade
Trademark	Sklar Edge™ TC

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • [Service@SklarCorp.com](mailto:Service@SklarCorp.com)