

610.756.7884

Home » Products » Sklar Edge TC Iris Scissors - 4-1/2"

## SKLAR EDGE TC IRIS SCISSORS - 4-1/2"

Product Number: 16-2505

4-1/2" Sklar Edge Tungsten Carbide Iris Scissors, Curved

Sklar Edge™ TC Iris Scissors provide an exceptionally durable and hard cutting edge. These scissors are designed to stay sharp for a significantly longer period of time than standard scissors. Sklar permanently bonds tungsten carbide inserts to the blades of the scissors which maximizes control and precision. Originally designed for ophthalmic procedures, Iris Scissors have become a multipurpose instrument used in various surgical environments. These fine, sharp scissors are ideal for detailed dissection of fine tissue and fine suture removal. In addition to ophthalmic procedures, Iris Scissors are now often used in emergency rooms and for OB/GYN and dermatologic procedures. Sklar Edge™ scissors can be identified by two gold coated finger rings. This product is curved with smooth, sharp/sharp tips and a length of 4-1/2 inches.



## **Additional Information**

Surname	Iris
Length	4-1/2"
Tip Configuration	Sharp/Sharp
Instrument Type	Iris Scissors
Curvature	Curved
Working Surface Style	Smooth
Handle	Finger Rings / Dual Gold
Material	Tungsten Carbide / Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade
Trademark	Sklar Edge™ TC

© 2016 Sklar Surgical Instruments. All Rights Reserved. 889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • Service@SklarCorp.com