



610.756.7884


[Home](#) » [Products](#) » [Sklar Edge TC Mayo-Stille Dissecting Scissors - 6-3/4"](#)

## SKLAR EDGE TC MAYO-STILLE DISSECTING SCISSORS - 6-3/4"

Product Number: 16-3557

6-3/4" Sklar Edge Tungsten Carbide Mayo-Stille Dissecting Scissors, Straight, Rounded Blades

Sklar Edge™ Tungsten Carbide Mayo-Stille Dissecting Scissors feature tungsten carbide inserts that are permanently bonded to the scissor blade. This provides an exceptionally hard and durable cutting edge and increases the overall life the instrument. Mayo-Stille Dissecting Scissors have blunt tips and are available curved or straight. The edges are rounded which permits deeper dissection of tough tissue. Sklar Edge™ scissors are differentiated from other scissors by the handle of two gold coated finger rings. This product is straight with smooth, blunt/blunt tips, and an overall length of 6-3/4 inches.



### Additional Information

Surname	Mayo-Stille
Length	6-3/4"
Tip Configuration	Blunt/Blunt
Instrument Type	Dissecting Scissors
Curvature	Straight
Working Surface Style	Smooth
Handle	Finger Rings / Dual Gold
Material	Tungsten Carbide / Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade
Trademark	Sklar Edge™ TC

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • [Service@SklarCorp.com](mailto:Service@SklarCorp.com)