



610.756.7884


[Home](#) » [Products](#) » [TC Mayo-Hegar Needle Holder](#)

TC MAYO-HEGAR NEEDLE HOLDER

Product Number: 21-8031

6" TC Mayo-Hegar Needle Holder Serrated Tungsten Carbide

The Sklar® TC Mayo-Hegar Needle Holder is a ratcheted, finger ring instrument available in a wide range of lengths, and used in multiple types of surgical procedures. Due to its versatility, it is one of the most commonly used surgical tools for both human and animal procedures. Shorter length products may be used in oral, gynecological, or dermatological procedures, as well as with small animals. Larger lengths of this product may be used for deep tissue suturing such as in bariatric procedures, or in large animal procedures. The Mayo-Hegar needle holder is commonly used with larger, heavier suture needles because of its wide jaws. Tungsten carbide inserts also maximize control and precision. This product is straight with cross-serrated jaws, and an overall length of 6 inches.



Additional Information

Surname	Mayo-Hegar
Length	6"
Tip Configuration	Blunt
Curvature	Straight
Working Surface Style	Cross-Serrated
Handle	Finger Rings / Dual Gold
Material	Tungsten Carbide / Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • Service@SklarCorp.com