



610.756.7884


[Home](#) » [Products](#) » [Spencer \(Littauer\) Stitch Scissors - 3-1/2"](#)

SPENCER (LITTAUER) STITCH SCISSORS - 3-1/2"

Product Number: 22-2835

3-1/2" Spencer Stitch Scissors, Smooth, Blunt/Blunt

Sklar® Spencer Stitch Scissors are primarily used for suture removal. These scissors have a small hook shaped tip on one blade that slides under sutures to slightly lift them before cutting for removal. This hook also holds the suture so it does not slip off the blade of the scissors prior to cutting. They are most commonly available in lengths of 3-1/2" and 4-1/2" and therefore typically have smaller, more delicate tips that are used for removing smaller sutures. They should not be used for cutting tissue. This product is straight with smooth, blunt/blunt tips and a length of 3-1/2 inches.



Additional Information

Surname	Spencer (Littauer)
Length	3-1/2"
Tip Configuration	Blunt/Blunt
Instrument Type	Scissors
Curvature	Straight
Working Surface Style	Smooth
Handle	Finger Rings
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • Service@SklarCorp.com