



610.756.7884


[Home](#) » [Products](#) » [Semken Tissue Forceps](#)

SEMKEN TISSUE FORCEPS

Product Number: 22-4354

6" Semken Tissue Forceps, Straight, 1x2 Teeth

Sklar® Semken Tissue Forceps are lightweight, narrow, toothed forceps used for handling primarily fine membranes or single layer tissues. This product is straight with 1x2 teeth, and a length of 6 inches.



Additional Information

| | |
|------------------------|-----------------------------------|
| Surname | Semken (sometimes spelled Semkin) |
| Length | 6" |
| Tip Configuration | 1x2 Teeth |
| Instrument Type | Tissue Forceps |
| Curvature | Straight |
| Pattern | Narrow |
| Handle | Thumb |
| Material | Stainless Steel |
| Disposable or Reusable | Reusable |
| Sterile or Non-Sterile | Non-Sterile |
| Latex or Latex-Free | Latex-Free |
| Grade | Premium OR-Grade |

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • Service@SklarCorp.com