



610.756.7884

[Home](#) » [Products](#) » [Sklarlite XD Dressing Forceps - 5-1/2"](#)

SKLARLITE XD DRESSING FORCEPS - 5-1/2"

Product Number: 23-2838

5-1/2" Sklarlite Extra Delicate Dressing Forceps

Sklarlite XD™ Dressing Forceps are thumb forceps used for holding dressing materials such as cotton and gauze during surgical procedures, changing dressings, or packing wounds. Sklarlite XD™ instruments are noticeably lighter and delicate than standard instruments which allow the surgeon the highest degree of precision with the least amount of fatigue. These instruments have been uniquely manufactured to be one of the lightest, most delicate general instruments available today. This product is straight with serrated tips and a length of 5-1/2 inches.



Additional Information

Length	5-1/2"
Instrument Type	Dressing Forceps
Curvature	Straight
Working Surface Style	Serrated
Pattern	Extra Light and Extra Delicate (SklarLite XD)
Handle	Thumb
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade
Trademark	Sklarlite XD™ Extra Light and Extra Delicate

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • Service@SklarCorp.com