



610.756.7884

[Home](#) » [Products](#) » [Probe with Eye - 5"](#)

PROBE WITH EYE - 5"

Product Number: 30-1950

5" Probe with Eye

Sklar's® Probe with Eye is a slender, slightly malleable, multipurpose instrument with a bulbous tip that can be used in a variety of procedures. It may be used to explore or gauge the depth of non-visible cavities, explore lesions, or manipulate tissue. This product is straight with an eye at the handle, and a length of 5 inches. The eye may be used for carrying ligatures.



Additional Information

Length	5"
Instrument Type	Probe
Size / Model	with Eye
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • Service@SklarCorp.com