

## 610.756.7863

Home » Products » Adson Tissue Forceps, 1x2 Teeth - 4-3/4"

## ADSON TISSUE FORCEPS, 1X2 TEETH - 4-3/4"

Product Number: 50-3050

4-3/4" Adson Tissue Forceps, 1x2 Teeth, Smooth Jaws, Straight, LiteGrip Fenestrated Handle

Sklar® LiteGrip™ Adson Tissue Forceps are for holding and manipulating delicate tissues. These forceps have a wide, flat thumb grasp area that is commonly serrated. The jaws are short and the tips are narrow. Sklar® LiteGrip™ instruments feature a lightweight, fenestrated handle. They are particularly useful during lengthy procedures when fatigue may arise. The finger holes on the handle also allow grip adjustments as needed by the surgeon. This product is straight with smooth jaws, 1X2 teeth, and a length of 4-3/4 inches.



## Additional Information

Surname	Adson
Length	4-3/4"
Tip Configuration	1x2 Teeth
Instrument Type	Tissue Forceps
Curvature	Straight
Working Surface Style	Smooth Jaws
Handle	Thumb
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade
Trademark	Sklar LiteGrip™ Fenestrated Handle

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7863 • Fax: 610.430.3941 • Service@SklarCorp.com