

610.756.7863

Home » Products » Adson TC Tissue Forceps, Tying Platform - 4-3/4"

ADSON TC TISSUE FORCEPS, TYING PLATFORM - 4-3/4"

Product Number: 51-3049

4-3/4" Adson Tungsten Carbide Tissue Forceps, Tying Platform, Serrated, 1x2 Teeth

Sklar® TC Adson Tissue Forceps are for holding and manipulating delicate tissues. Sklar permanently bonds tungsten carbide inserts to the instrument, which provides maximum control and precision for the user, as well as increasing the length of time the instrument will remain sharp. These forceps have a wide, flat thumb grasp area that is commonly slightly grooved for an improved grip. The jaws are short and the tips are narrow. These TC Adson tissue forceps, identified by gold handles, are straight with cross-serrated tips, 1 X 2 teeth, and a length of 4-3/4 inches.



Additional Information

Surname	Adson
Length	4-3/4"
Tip Configuration	1x2 Teeth
Instrument Type	Tissue Forceps
Curvature	Straight
Working Surface Style	Cross-Serrated Tying Platform
Handle	Thumb
Material	Tungsten Carbide / Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved. 889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7863 • Fax: 610.430.3941 • Service@SklarCorp.com