

610.756.7863

Home » Products » Senn Retractor

SENN RETRACTOR

Product Number: 60-6162

6-1/4" Senn Retractor, Double End, 1-3 Prong, Sharp, 1-Flat Blade, Blunt

Sklar's® Senn Retractor is a handheld, double-ended retractor used to retract primarily surface tissue. It is often used in plastic surgery, small bone and joint procedures, or thyroidectomy and dissection of neck tissue. This instrument has one angled, blunt end and, at the opposite end, facing the opposite direction has a three prong, rake tip. The prongs may be sharp or blunt. This product has a straight, flat handle with one squared off blunt end and a second rake end with sharp tips. The overall length is 6-1/4 inches.



Additional Information

Surname	Senn
Length	6-1/4"
Tip Configuration	3 Prongs Sharp
Instrument Type	Retractor
Curvature	Straight
Working Surface Style	Smooth
Handle	Smooth Flat
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved. 889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7863 • Fax: 610.430.3941 • Service@SklarCorp.com