



610.756.7863

[Home](#) » [Products](#) » [Hegar Uterine Dilator](#)

## HEGAR UTERINE DILATOR

Product Number: 90-4804

Hegar Uterine Dilator Double End 3mm/4mm

Sklar's® Hegar Double-Ended Uterine Dilator is a manual, probe-style instrument used to dilate the cervix. They are solid, nonflexible, and available in multiple diameters. These dilators are slightly curved in opposite directions from the center and have smooth, rounded tips. This product is the 3mm/4mm dilator. The overall length is 8 inches.



### Additional Information

Surname	Hegar
Instrument Type	Uterine Dilator
Curvature	Curved
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7863 • Fax: 610.430.3941 • [Service@SklarCorp.com](mailto:Service@SklarCorp.com)