



610.756.7865


[Home](#) » [Products](#) » [Surgi-OR Rochester-Pean Forceps](#)

SURGI-OR ROCHESTER-PEAN FORCEPS

Product Number: 95-470

5-1/2" Surgi-OR Rochester-Pean Forceps Straight

Surgi-OR™ Rochester-Pean Hemostatic Forceps are used for clamping larger tissue and vessels for hemostasis. The full horizontal serrations and availability in multiple lengths make these forceps a versatile instrument used in multiple procedures. They are available curved or straight. Hand-crafted in Pakistan, Surgi-OR™ mid-grade instruments are forged from the highest quality Japanese steel and provide reliable instrumentation for the operating room and other surgical centers. Surgi-OR™ instruments are a cost-effective line of instruments that must meet Sklar's quality standards. This product is straight with full, horizontally serrated jaws, and a length of 5-1/2 inches.



Additional Information

Surname	Rochester-Pean
Length	5-1/2"
Instrument Type	Hemostatic Forceps
Curvature	Straight
Working Surface Style	Serrated
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Mid-Grade (Pakistan)

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7865 • Fax: 610.430.3941 • Service@SklarCorp.com