



610.756.7865

[Home](#) » [Products](#) » [Surgi-OR Plain Splinter Forceps - 4-1/2"](#)

## SURGI-OR PLAIN SPLINTER FORCEPS - 4-1/2"

Product Number: 95-781

4-1/2" Surgi-OR Plain Splinter Forceps, Straight, Serrated

Surgi-OR™ Plain Splinter Forceps are multi-purpose thumb forceps typically used for grasping and removing small foreign objects with their fine, pointed tips. They may also be used for holding bandages, sutures, or other items. Hand-crafted in Pakistan, Surgi-OR™ mid-grade instruments are forged from the highest quality Japanese steel and provide reliable instrumentation for the operating room and other surgical centers. Surgi-OR™ instruments are a cost-effective line of instruments that must meet Sklar's quality standards. This product is straight with serrated jaws and pointed tips. It has a length of 4-1/2 inches.



### Additional Information

Length	4-1/2"
Instrument Type	Splinter Forceps
Curvature	Straight
Working Surface Style	Serrated
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Mid-Grade (Pakistan)

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7865 • Fax: 610.430.3941 • [Service@SklarCorp.com](mailto:Service@SklarCorp.com)