



610.756.7865


[Home](#) » [Products](#) » [Nail Nipper - Double Spring with Lock, 5-1/2"](#)

## NAIL NIPPER - DOUBLE SPRING WITH LOCK, 5-1/2"

Product Number: 97-1056

5-1/2" Nail Nipper, Straight Back, Double Spring with Lock

Sklar® Nail Nippers are used for removing portions of the toe nail or the nail plate. This product features double spring, grooved, plier handles with a locking mechanism that is released upon squeezing the handles. It is 5-1/2 inches long and features concave blades with a straight back.



### Additional Information

Length	5-1/2"
Instrument Type	Nail Nipper
Curvature	Straight
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7865 • Fax: 610.430.3941 • [Service@SklarCorp.com](mailto:Service@SklarCorp.com)