



610.756.7865


[Home](#) » [Products](#) » [Merit Operating Scissors - B/B, 5-1/2"](#)

## MERIT OPERATING SCISSORS - B/B, 5-1/2"

Product Number: 97-278

5-1/2" Merit Operating Scissors, Straight, Blunt/Blunt

Merit™ Operating Scissors, also known as Surgical Scissors, are a versatile surgical instrument. With three available tip configurations (sharp/sharp, blunt/blunt, or sharp/blunt), varying curvatures and multiple lengths, these scissors are used for general, multipurpose cutting and dissecting. Sharp Operating Scissors cut tougher tissue and are frequently used for dissecting, whereas blunt Operating Scissors are more often used for cutting materials such as sutures, tubing, or gauze. Merit™ instruments are Sklar's office-grade line of products and most commonly utilized for routine and in-office minor procedures. These cost effective, high quality instruments are crafted from Pakistan steel and come with a 10 year guarantee. Merit™ instruments deliver an economical solution for physicians and practices that require quality, durable instruments. This product is straight with smooth, blunt/blunt tips and a length of 5-1/2 inches.



### Additional Information

Length	5-1/2"
Tip Configuration	Blunt/Blunt
Instrument Type	Operating Scissors
Curvature	Straight
Working Surface Style	Smooth
Handle	Finger Rings
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Physician Grade (Pakistan)

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7865 • Fax: 610.430.3941 • [Service@SklarCorp.com](mailto:Service@SklarCorp.com)