



610.756.7865


[Home](#) » [Products](#) » [Merit Strabismus Scissors](#)

## MERIT STRABISMUS SCISSORS

Product Number: 97-361

4-1/2" Merit Strabismus Scissors Curved

Merit™ instruments are Sklar's office-grade line of products and most commonly utilized for routine and in-office minor procedures. These cost effective, high quality instruments are crafted from Pakistan steel and come with a 10 year guarantee. Merit™ instruments deliver an economical solution for physicians and practices that require quality, durable instruments. This product is straight with serrated tips and a length of 5-1/2 inches.



### Additional Information

Length	4-1/2"
Tip Configuration	Blunt/Blunt
Instrument Type	Scissors
Curvature	Curved
Working Surface Style	Smooth
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Physician Grade (Pakistan)

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7865 • Fax: 610.430.3941 • [Service@SklarCorp.com](mailto:Service@SklarCorp.com)