



610.756.7879

[Home](#) » [Products](#) » [Merit Adson Dressing Forceps - 4-3/4"](#)

## MERIT ADSON DRESSING FORCEPS - 4-3/4"

Product Number: 97-771

4-3/4" Merit Adson Dressing Forceps, Serrated, Straight



### Additional Information

Surname	Adson
Length	4-3/4"
Instrument Type	Dressing Forceps
Curvature	Straight
Working Surface Style	Serrated
Handle	Thumb
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Physician Grade (Pakistan)

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7879 • Fax: 610.430.3941 • [Service@SklarCorp.com](mailto:Service@SklarCorp.com)