

## 610.756.7879

Home » Products » Merit Plain Splinter Forceps - 3-1/2"

## MERIT PLAIN SPLINTER FORCEPS - 3-1/2"

Product Number: 97-777

3-1/2" Merit Plain Splinter Forceps, Straight, Serrated

Merit ™ Plain Splinter Forceps are multi-purpose thumb forceps typically used for grasping and removing small foreign objects with their fine, pointed tips. They may also be used for holding bandages, sutures, or other items. Merit™ instruments are Sklar's office-grade line of products and most commonly utilized for routine and in-office minor procedures. These cost effective, high quality instruments are crafted from Pakistan steel and come with a 10 year guarantee. Merit™ instruments deliver an economical solution for physicians and practices that require quality, durable instruments. This product is straight with serrated jaws and pointed tips. It has a length of 3-1/2 inches.



## Additional Information

Length	3-1/2"
Instrument Type	Splinter Forceps
Curvature	Straight
Working Surface Style	Serrated
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Physician Grade (Pakistan)

© 2016 Sklar Surgical Instruments. All Rights Reserved. 889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7879 • Fax: 610.430.3941 • Service@SklarCorp.com