



610.756.7879


[Home](#) » [Products](#) » [Merit Carmalt Splinter Forceps - 4-1/4"](#)

## MERIT CARMALT SPLINTER FORCEPS - 4-1/4"

Product Number: 97-791

4-1/4" Merit Carmalt Splinter Forceps, Straight, Fine Point

Merit™ Carmalt Splinter Forceps are thumb forceps commonly used in dermatological procedures. They feature an arrowhead style tip with a central groove. They may also be known as Hunter forceps. Merit™ instruments are Sklar's office-grade line of products and most commonly utilized for routine and in-office minor procedures. This product is straight with serrated, fine point tips, and a length of 4-1/4 inches.



### Additional Information

Surname	Carmalt
Length	4-1/4"
Tip Configuration	Fine Point
Instrument Type	Splinter Forceps
Curvature	Straight
Working Surface Style	Serrated
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Physician Grade (Pakistan)

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7879 • Fax: 610.430.3941 • [Service@SklarCorp.com](mailto:Service@SklarCorp.com)