



610.756.7879


[Home](#) » [Products](#) » [Merit Iris Scissors](#)

## MERIT IRIS SCISSORS

Product Number: 98-109

4-1/2" Merit Iris Scissors, Curved, Smooth, Sharp/Sharp

Merit™ Iris Scissors were originally designed for ophthalmic procedures, but have become a multipurpose instrument used in various surgical environments. These fine, sharp scissors are ideal for detailed dissection of fine tissue and fine suture removal. In addition to ophthalmic procedures, Iris Scissors are now often used for multiple procedures. Merit™ instruments are Sklar's office-grade line of products and most commonly utilized for routine and in-office minor procedures. This product is curved with smooth, sharp/sharp tips and a length of 4-1/2 inches.



### Additional Information

Surname	Iris
Length	4-1/2"
Tip Configuration	Sharp/Sharp
Instrument Type	Iris Scissors
Curvature	Curved
Working Surface Style	Smooth
Handle	Finger Rings
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Physician Grade (Pakistan)

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7879 • Fax: 610.430.3941 • [Service@SklarCorp.com](mailto:Service@SklarCorp.com)