



610.756.7884

[Home](#) » [Products](#) » [Cushing Brain Forceps](#)

CUSHING BRAIN FORCEPS

Product Number: 22-6770

7" Cushing Brain Forceps Delicate Serrated



Additional Information

Surname	Cushing
Length	7"
Instrument Type	Tissue Forceps
Curvature	Straight
Working Surface Style	Serrated
Pattern	Delicate
Handle	Thumb
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • Service@SklarCorp.com