



610.756.7879


[Home](#) » [Products](#) » [Merit Micro Alligator Forceps](#)

MERIT MICRO ALLIGATOR FORCEPS

Product Number: 98-239

2-3/4" Merit Micro Alligator Forceps Serrated

Merit™ Micro Alligator Ear Forceps are non-ratcheted, finger ring forceps used for manipulating and removing tissue from the ear canal or middle ear. They may also be used for packing during grafting procedures or as nasal forceps. Merit™ instruments are Sklar's office-grade line of products and most commonly utilized for routine and in-office minor procedures. These cost effective, high quality instruments are crafted from Pakistan steel and come with a 10 year guarantee. Merit™ instruments deliver an economical solution for physicians and practices that require quality, durable instruments. This product has 2-3/4 inch straight shanks. At the box lock the shaft extends 3-1/4 inches to 6mm serrated jaws.



Additional Information

Surname	Alligator
Length	2-3/4"
Instrument Type	Ear/Nasal Forceps
Curvature	Straight
Working Surface Style	Serrated
Pattern	Very Delicate
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Physician Grade (Pakistan)

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7879 • Fax: 610.430.3941 • Service@SklarCorp.com