



610.756.7879

[Home](#) » [Products](#) » [Merit Vienna Nasal Speculum - Adult](#)

MERIT VIENNA NASAL SPECULUM - ADULT

Product Number: 98-258

Merit Vienna Nasal Speculum, Adult

The Merit™ Vienna Nasal Speculum is an intranasal instrument used to retract and expose nasal tissue for examination. Merit™ instruments are Sklar's office-grade line of products and most commonly utilized for routine and in-office minor procedures. These cost effective, high quality instruments are crafted from Pakistan steel and come with a 10 year guarantee. Merit™ instruments deliver an economical solution for physicians and practices that require quality, durable instruments. This product is the adult size speculum.



Additional Information

Surname	Vienna (also known as Wien Model, Mod. Wien)
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Physician Grade (Pakistan)

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7879 • Fax: 610.430.3941 • Service@SklarCorp.com