

Keyes Ring Cutter | Sklar Surgical Instruments



Home » Products » Keyes Ring Cutter

KEYES RING CUTTER Product Number: 06-4150

Keyes Ring Cutter Stainless Steel. For use on soft metals only. Use Blades 06-4151.

Sklar's® Keyes Ring Cutter is a speciality instrument specifically used for cutting and removing rings from a traumatized or otherwise inflamed finger. It features a thumb screw that turns a sharp, toothed, rotary blade that cuts the ring. The plier handle operates a lower jaw that slides under the ring and protects the finger from the cutting while holding the blade onto the surface of the ring. The Keyes ring cutter is best used on soft metals.



Additional Information

Surname	Keyes
Instrument Type	Ring Cutter
Working Surface Style	Sharp Cutting Edge
Handle	Pliers
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7888 • Fax: 610.430.3941 • Service@SklarCorp.com