



610.756.7884


[Home](#) » [Products](#) » [Stille Plaster Shears](#)

STILLE PLASTER SHEARS

Product Number: 11-1615

15" Stille Plaster Shears Two Handed

Sklar's® Stille Plaster Shears are large 2 handed shears used to remove heavy or large plaster casts. They have round handles to facilitate ease in the cast removal.



Additional Information

Surname	Stille
Length	15"
Instrument Type	Plaster Shears
Curvature	Offset
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • Service@SklarCorp.com