



610.756.7884

[Home](#) » [Products](#) » [Hennig Plaster Spreader](#)

HENNIG PLASTER SPREADER

Product Number: 11-1910

11" Hennig Plaster Spreader

Sklar's® Hennig Plaster Spreader is used to aid in the removal of plaster casts. It has a large plier shape handle and two large, flat, rectangular blades with horizontal serrations on the outside of the blade. Once a cut is made in the plaster cast, the blades of the spreader are closed by opening the handle and are inserted to a section of the cast. By squeezing the handle, the blades spread open, there by spreading apart sections of the cast. After multiple sections of the cast are spread open, the cast may be removed. The blades are approximately 2-1/8 inches wide and the overall length of this instrument is 11 inches.



Additional Information

Surname	Hennig
Length	11"
Tip Configuration	2 Flat Blades Outside Serrated
Instrument Type	Plaster Spreader
Curvature	Straight
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • Service@SklarCorp.com