



610.756.7884


[Home](#) » [Products](#) » [Sklar Edge TC Metzenbaum-Lahey Dissecting Scissors - 6"](#)

SKLAR EDGE TC METZENBAUM-LAHEY DISSECTING SCISSORS - 6"

Product Number: 16-1910

6" Sklar Edge Tungsten Carbide Metzenbaum-Lahey Dissecting Scissors, Curved

Sklar Edge™ Tungsten Carbide Metzenbaum-Lahey Dissecting Scissors are most commonly used for dissecting and cutting tissue. Metzenbaum-Lahey scissors tend to be shorter than traditional Metzenbaum scissors. Sklar Edge™ Tungsten Carbide Scissors provide an exceptionally hard and durable cutting edge. These scissors are designed to stay sharp for a significantly longer period of time than other non tungsten carbide scissors. Sklar permanently bonds the tungsten carbide inserts to the blade of the scissors which maximizes control and precision as well. Sklar Edge™ scissors are differentiated from other scissors by the handle of two gold coated finger rings. This product has curved, smooth, blunt tips, and an overall length of 6 inches.



Additional Information

Surname	Metzenbaum-Lahey
Length	6"
Tip Configuration	Blunt/Blunt
Instrument Type	Dissecting Scissors
Curvature	Curved
Working Surface Style	Smooth
Handle	Finger Rings / Dual Gold
Material	Tungsten Carbide / Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade
Trademark	Sklar Edge™ TC

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • Service@SklarCorp.com