



610.756.7884

[Home](#) » [Products](#) » [Sklar Edge TC Metzenbaum Dissecting Scissors - 7"](#)

SKLAR EDGE TC METZENBAUM DISSECTING SCISSORS - 7"

Product Number: 16-2010

7" Sklar Edge Tungsten Carbide Metzenbaum Dissecting Scissors, Straight

Sklar Edge™ Tungsten Carbide Metzenbaum Dissecting Scissors are most commonly used for dissecting and cutting tissue. Sklar Edge™ Tungsten Carbide Scissors provide an exceptionally durable and hard cutting edge. These scissors are designed to stay sharp for a significantly longer period of time than standard scissors. Sklar permanently bonds tungsten carbide inserts to the blades of the scissors which maximizes control and precision. Sklar Edge™ scissors can be identified by two gold coated finger rings. This product is straight with smooth edges, blunt tips, and an overall length of 7 inches.



Additional Information

Surname	Metzenbaum
Length	7"
Tip Configuration	Blunt/Blunt
Instrument Type	Dissecting Scissors
Curvature	Straight
Working Surface Style	Smooth
Handle	Finger Rings / Dual Gold
Material	Tungsten Carbide / Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade
Trademark	Sklar Edge™ TC

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • Service@SklarCorp.com