



610.756.7884

[Home](#) » [Products](#) » [Sklar Edge TC Metzenbaum-Nelson Bariatric](#)

[Dissecting Scissors - 9"](#)

SKLAR EDGE TC METZENBAUM-NELSON BARIATRIC DISSECTING SCISSORS - 9"

Product Number: 16-2310

9" Sklar Edge Tungsten Carbide Metzenbaum-Nelson Bariatric Dissecting Scissors, Curved

Sklar Edge™ Tungsten Carbide Metzenbaum-Nelson Bariatric Scissors are long dissecting scissors most commonly used in Bariatric Surgeries such as gastric bypass, gastric banding, and sleeve gastrectomy. They are also used in Nephrectomy, Cystectomy, and Prostatectomy procedures to dissect and cut tissue. Metzenbaum-Nelson scissors tend to be longer than standard Metzenbaum scissors which make them ideal for use in thoracic procedures as well. These scissors are not recommended for cutting sutures or drains. Sklar Edge™ TC Scissors provide an exceptionally durable and hard cutting edge. These scissors are designed to stay sharp for a significantly longer period of time than standard scissors. Sklar permanently bonds tungsten carbide inserts to the blades of the scissors which maximizes control and precision. Sklar Edge™ scissors can be identified by two gold coated finger rings. This product has curved, smooth blades with blunt tips, and an overall length of 9 inches.



Additional Information

Surname	Metzenbaum-Nelson
Length	9"
Tip Configuration	Blunt/Blunt
Instrument Type	Dissecting Scissors
Curvature	Curved
Working Surface Style	Smooth
Handle	Finger Rings / Dual Gold
Material	Tungsten Carbide / Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade
Trademark	Sklar Edge™ TC

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • Service@SklarCorp.com