



610.756.7884

[Home](#) » [Products](#) » [Rochester-Pean Forceps - 8"](#)

ROCHESTER-PEAN FORCEPS - 8"

Product Number: 17-2480

8" Rochester-Pean Forceps Straight

Sklar® Rochester-Pean Hemostatic Forceps are used for clamping larger tissue and vessels for hemostasis. The full horizontal serrations and availability in multiple lengths make these forceps a versatile instrument used in multiple procedures. They are available curved or straight. This product is straight with full, horizontally serrated jaws, and a length of 8 inches.



Additional Information

Surname	Rochester-Pean
Length	8"
Instrument Type	Hemostatic Forceps
Curvature	Straight
Working Surface Style	Serrated
Handle	Finger Rings
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • Service@SklarCorp.com