



610.756.7884

[Home](#) » [Products](#) » [Standard Fine Splinter Forceps - 4"](#)

STANDARD FINE SPLINTER FORCEPS - 4"

Product Number: 19-3040

4" Standard Fine Splinter Forceps

Sklar® Plain Splinter Forceps are multi-purpose thumb forceps typically used for grasping and removing small foreign objects with their fine, pointed tips. They may also be used for holding bandages, sutures, or other items. This product is straight with serrated jaws and fine, pointed tips. It has a length of 4 inches.



Additional Information

| | |
|------------------------|------------------|
| Length | 4" |
| Tip Configuration | Sharp |
| Instrument Type | Splinter Forceps |
| Curvature | Straight |
| Working Surface Style | Serrated |
| Pattern | Fine |
| Handle | Thumb |
| Material | Stainless Steel |
| Disposable or Reusable | Reusable |
| Sterile or Non-Sterile | Non-Sterile |
| Latex or Latex-Free | Latex-Free |
| Grade | Premium OR-Grade |

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • Service@SklarCorp.com