

610.756.7884

Home » Products » Olsen-Hegar Needle Holder

OLSEN-HEGAR NEEDLE HOLDER

Product Number: 20-2755

5-1/2" Olsen-Hegar Needle Holder, Cross-Serrated with Groove

Sklar's® Olsen-Hegar Needle Holder is a ratcheted, finger ring instrument that also has built in scissors. It is used with intermediate weight sutures, commonly in plastic surgery procedures. This product is straight with cross-serrated tips, and a central groove. It is 5-1/2 inches long.



Additional Information

Surname	Olsen-Hegar
Length	5-1/2"
Tip Configuration	Blunt
Instrument Type	Needle Holder
Curvature	Straight
Working Surface Style	Cross-Serrated with Groove
Handle	Finger Rings
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved. 889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • Service@SklarCorp.com