



610.756.7884

[Home](#) » [Products](#) » [TC Olsen-Hegar Needle Holder](#)

## TC OLSEN-HEGAR NEEDLE HOLDER

Product Number: 21-8027

5-1/2" TC Olsen-Hegar Needle Holder Serrated Tungsten Carbide

Sklar's® Olsen-Hegar TC Needle Holder is a ratcheted, finger ring instrument that also has built in scissors. It is used with intermediate weight sutures, commonly in plastic surgery procedures. Sklar permanently bonds tungsten carbide inserts to the instrument, which maximizes control and precision. This product is straight with cross-serrated tips, and an overall length of 5-1/2 inches.



### Additional Information

Surname	Olsen-Hegar
Length	5-1/2"
Tip Configuration	Blunt
Curvature	Straight
Working Surface Style	Cross-Serrated
Handle	Finger Rings / Dual Gold
Material	Tungsten Carbide / Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • [Service@SklarCorp.com](mailto:Service@SklarCorp.com)