



610.756.7884

[Home](#) » [Products](#) » [Allis Tissue Forceps - 6"](#)

ALLIS TISSUE FORCEPS - 6"

Product Number: 36-2060

6" Allis Tissue Forceps, 3x4 Teeth

Sklar® Allis Tissue Forceps are finger ring, ratcheted forceps used for grasping organs and slippery or dense tissue during surgical procedures. They are commonly used with tonsil, vaginal, breast, and thyroid tissue. Available straight or angled, this instrument is offered in a variety of lengths. The tip end style varies to meet the needs of the procedure or preferences of the surgeon. This product is straight with 3x4 teeth and an overall length of 6 inches.



Additional Information

Surname	Allis
Length	6"
Tip Configuration	3x4 Teeth
Instrument Type	Tissue Grasping Forceps
Curvature	Straight
Handle	Finger Rings
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • Service@SklarCorp.com