



610.756.7884

[Home](#) » [Products](#) » [Guthrie Fixation Hook - 2 Large Prongs, 5"](#)

GUTHRIE FIXATION HOOK - 2 LARGE PRONGS, 5"

Product Number: 47-2650

5" Guthrie Fixation Hook, 2 Large Prongs, Sharp

Sklar's® Guthrie Hook is a double-prong tissue hook used for handling and retracting tissue. They are often used in dermatologic and plastic surgery procedures but the tiny hooks also make this retractor useful in ophthalmologic and intranasal procedures. It is straight with two sharp, hooked, prong tips. This product is considered the large hook with the prongs 3.3mm to 3.5mm apart and an overall length of 5 inches.



Additional Information

Surname	Guthrie
Length	5"
Tip Configuration	2 Prongs Sharp
Instrument Type	Hook
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • Service@SklarCorp.com