



610.756.7884


[Home](#) » [Products](#) » [Adson Tissue Forceps - Serrated, 4-3/4"](#)

ADSON TISSUE FORCEPS - SERRATED, 4-3/4"

Product Number: 50-3048

4-3/4" Adson Tissue Forceps, Cross-Serrated Jaws, 1x2 Teeth, Straight

Sklar® Adson Tissue Forceps are for holding and manipulating delicate tissues. These forceps have a wide, flat thumb grasp area that is commonly serrated. The jaws are short and the tips are narrow. This product is straight with cross-serrated jaws, 1x2 teeth, and a length of 4-3/4 inches.



Additional Information

Surname	Adson
Length	4-3/4"
Tip Configuration	1x2 Teeth
Instrument Type	Tissue Forceps
Curvature	Straight
Working Surface Style	Cross-Serrated Jaws
Handle	Thumb
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7884 • Fax: 610.430.3941 • Service@SklarCorp.com