



610.756.7863

[Home](#) » [Products](#) » [Adson Tissue Forceps, Tying Platform, 1x2 Teeth](#)


- 4-3/4"

## ADSON TISSUE FORCEPS, TYING PLATFORM, 1X2 TEETH - 4-3/4"

Product Number: 50-3049

4-3/4" Adson Tissue Forceps, Tying Platform, 1x2 Teeth, Straight

Sklar® Adson Tissue Forceps are for holding and manipulating delicate tissues. These forceps have a wide, flat thumb grasp area that is commonly serrated. The jaws are short and the tips are narrow. This product is straight with smooth jaws, 1x2 teeth, and a length of 4-3/4 inches.



### Additional Information

Surname	Adson
Length	4-3/4"
Tip Configuration	1x2 Teeth
Instrument Type	Tissue Forceps
Curvature	Straight
Working Surface Style	Smooth Tying Platform
Handle	Thumb
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7863 • Fax: 610.430.3941 • [Service@SklarCorp.com](mailto:Service@SklarCorp.com)