



610.756.7863

[Home](#) » [Products](#) » [Adson Dressing Forceps, Serrated - 4-3/4"](#)

## ADSON DRESSING FORCEPS, SERRATED - 4-3/4"

Product Number: 50-3148

4-3/4" Adson Dressing Forceps, Serrated, LiteGrip Fenestrated Handle

Sklar® Adson Dressing Forceps are thumb forceps used for holding dressing materials such as cotton and gauze during surgical procedures, changing dressings, or packing wounds. They have a wide thumb grasp for increased precision and control. Sklar® LiteGrip™ instruments feature a lightweight, fenestrated handle. They are particularly useful during lengthy procedures when fatigue may arise. The finger holes on the handle also allow grip adjustments as needed by the surgeon. This product is straight with serrated tips and a length of 4-3/4 inches.



### Additional Information

Surname	Adson
Length	4-3/4"
Instrument Type	Dressing Forceps
Curvature	Straight
Working Surface Style	Serrated
Handle	Thumb
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade
Trademark	Sklar LiteGrip™ Fenestrated Handle

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7863 • Fax: 610.430.3941 • [Service@SklarCorp.com](mailto:Service@SklarCorp.com)