



610.756.7863

[Home](#) » [Products](#) » [Bergh Cilia Forceps](#)

BERGH CILIA FORCEPS

Product Number: 66-1635

3-1/2" Bergh Cilia Forceps

Sklar® Bergh Cilia Forceps are heavy spring thumb forceps used during ophthalmologic procedures involving the eyelid, commonly, blepharitis marginalis. These forceps are horizontally serrated and used for grasping and removing cilia, or eyelashes. The flat, grooved handle provides a secure grip for the user. This product has horizontal serrations and an overall length of 3-1/2 inches.



Additional Information

Surname	Bergh
Length	3-1/2"
Instrument Type	Tissue Forceps
Curvature	Straight
Working Surface Style	Horizontally Serrated
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7863 • Fax: 610.430.3941 • Service@SklarCorp.com