



610.756.7863


[Home](#) » [Products](#) » [Hartmann Alligator Ear Forceps - Serrated, 3-1/2"](#)

HARTMANN ALLIGATOR EAR FORCEPS - SERRATED, 3-1/2"

Product Number: 67-1035

3-1/2" Hartmann Alligator Ear Forceps, Horizontally Serrated

Sklar® Hartmann Alligator Forceps are an ENT instrument commonly used in myringotomy procedures to insert or remove tubes or dressings and packing. The serrated tips securely hold and remove foreign objects, granulation tissue, and cerumen. This product has angled shanks, straight jaws, horizontally serrated tips and a length of 3-1/2 inches.



Additional Information

Surname	Hartmann / Alligator
Length	3-1/2"
Instrument Type	Forceps
Curvature	Straight Jaws with Angled Shanks
Working Surface Style	Horizontally Serrated
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7863 • Fax: 610.430.3941 • Service@SklarCorp.com