



610.756.7863


[Home](#) » [Products](#) » [Micro Alligator Ear Forceps](#)

MICRO ALLIGATOR EAR FORCEPS

Product Number: 67-8006

3-1/4" Micro Alligator Ear Forceps Delicate Straight 4mm Jaw

Sklar® Alligator Ear Forceps are non-ratcheted, finger ring ENT forceps primarily used for manipulating and removing tissue from the ear canal or middle ear. They may also be used for packing during grafting procedures or as nasal forceps. This product has 2-3/4 inch straight shanks which then angle to a 3-1/4 inch shaft. The jaws are 4mm long and serrated.



Additional Information

Surname	Alligator
Length	3-1/4"
Tip Dimensions	4mm
Instrument Type	Ear/Nasal Forceps
Curvature	Straight
Working Surface Style	Serrated
Pattern	Delicate
Handle	Finger Rings
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7863 • Fax: 610.430.3941 • Service@SklarCorp.com