



610.756.7863

[Home](#) » [Products](#) » [Wilson Sutureless Vasectomy Set](#)

WILSON SUTURELESS VASECTOMY SET

Product Number: 85-6220

Wilson Sutureless Vasectomy Set of 2, Ring Tip. Consists of 5-3/4" Sharply Pointed Forceps to Puncture the Scrotal Skin and a 5-1/2" Elevation Clamp to Encircle and Elevate the Vas. - Set of 2

Sklar's® Wilson Vasectomy Set consists of 2 instruments used to complete a non-incision, sutureless vasectomy. It includes 5-1/2 inch Curved Vas Piercing Forceps with sharp tips used to pierce the scrotal skin; and a 5-3/4 inch Vasectomy Clamp with ring tips to encircle and elevate the Vas Deferens. Both of the instruments have finger ring handles and are ratcheted. The tubes are then tied off, cauterized, or blocked.



Additional Information

Surname	Wilson
Instrument Type	Procedure Kit
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7863 • Fax: 610.430.3941 • Service@SklarCorp.com