



610.756.7863

[Home](#) » [Products](#) » [Vasectomy Fixation Clamp - Ring Tip, 5-3/4"](#)

VASECTOMY FIXATION CLAMP - RING TIP, 5-3/4"

Product Number: 85-6221

No-Scalpel Vasectomy Fixator Clamp Forceps, Round Ring Tip, 5-3/4" Overall length

Sklar's® Vasectomy Clamp with ring tips is a ratcheted, finger ring instrument used to encircle and elevate the Vas Deferens in non-suture vasectomy procedures. It is used in conjunction with vasectomy piercing forceps. This product is straight with ring tips and 5-3/4 inches long.



Additional Information

Length	5-3/4"
Tip Configuration	Ring Tip
Instrument Type	Hemostatic Forceps
Curvature	Straight
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7863 • Fax: 610.430.3941 • Service@SklarCorp.com