

## 610.756.7863

Home » Products » Foerster Sponge Forceps - Straight, Serrated, 9-1/2"

## FOERSTER SPONGE FORCEPS - STRAIGHT, SERRATED, 9-1/2"

Product Number: 87-2195

9-1/2" Foerster Sponge Forceps, Serrated, Straight

Sklar® Foerster Sponge Forceps, which may also be known as Foerster-Ballenger Sponge Forceps, are commonly used in surgical procedures to hold gauze squares that will absorb excess fluids and blood from the surgical site or hemostasis. The ratcheted locking mechanism of these finger ring forceps securely grip the gauze and hold it in place. Foerster sponge forceps are often used in bariatric surgeries because of their longer length, and the curved version may be used for removing the placenta in OB/GYN procedures. This product is straight with oval, fenestrated, serrated tips and a length of 9-1/2 inches.



## **Additional Information**

Surname	Foerster-Ballenger
Length	9-1/2"
Tip Configuration	Fenestrated
Instrument Type	Sponge Forceps
Curvature	Straight
Working Surface Style	Serrated
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved. 889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7863 • Fax: 610.430.3941 • Service@SklarCorp.com