



610.756.7863

[Home](#) » [Products](#) » [Simpson Uterine Sound](#)

SIMPSON UTERINE SOUND

Product Number: 90-5413

12-1/2" Simpson Uterine Sound Centimeter Graduation Silver Plated

Sklar's® Simpson Uterine Sound is a malleable, silver plated OB/GYN instrument used for probing and dilating the uterus, as well as measuring the length of the cervical canal. It is round and thin with a small nodule at the tip and a thumb grip handle. It is very similar to the Sims uterine sound. This product has a gently curved tip, graduated centimeter markings, and an overall length of 12-1/2 inches.



Additional Information

Surname	Simpson
Length	12-1/2"
Instrument Type	Uterine Sound
Curvature	Curved
Material	Silver Plated
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7863 • Fax: 610.430.3941 • Service@SklarCorp.com