



610.756.7863


[Home](#) » [Products](#) » [Teale Vulsellum Forceps](#)

TEALE VULSELLUM FORCEPS

Product Number: 91-2695

9" Teale Vulsellum Forceps Curved Sideways 3x4 Teeth

Sklar's® Teale Vulsellum Forceps is a ratcheted, finger ring, OB/GYN instrument used for gripping the cervix or uterus during procedures such as episiotomy, hysterectomy, or dilation & curettage. These forceps have 3x4 teeth that provide a secure grip with less trauma than single toothed tenaculum forceps, which bite more deeply. Teale vulsellum forceps have long shanks and long jaws, and are available straight, angled, or curved sideways. This product is curved sideways and has an overall length of 9 inches.



Additional Information

Surname	Teale
Length	9"
Tip Configuration	3x4 Teeth
Instrument Type	Grasping Forceps
Curvature	Curved (Sideways)
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7863 • Fax: 610.430.3941 • Service@SklarCorp.com