



610.756.7865


[Home](#) » [Products](#) » [Surgi-OR Strabismus Scissors](#)

SURGI-OR STRABISMUS SCISSORS

Product Number: 95-361

4-1/2" Surgi-OR Strabismus Scissors Curved

Surgi-OR™ Strabismus Scissors are an ophthalmic dissecting instrument primarily used during procedures for the treatment of Strabismus, otherwise known as "cross-eyes." Hand-crafted in Pakistan, these mid-grade instruments are forged from the highest quality Japanese steel and provide reliable instrumentation for the operating room and other surgical centers. Surgi-OR™ instruments are a cost-effective line of instruments that must meet Sklar's quality standards. This product is curved with smooth, blunt/blunt tips and a length of 4-1/2 inches.



Additional Information

Length	4-1/2"
Tip Configuration	Blunt/Blunt
Instrument Type	Scissors
Curvature	Curved
Working Surface Style	Smooth
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Mid-Grade (Pakistan)

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7865 • Fax: 610.430.3941 • Service@SklarCorp.com