



610.756.7865


[Home](#) » [Products](#) » [Ingrown Nail Splitting Scissors](#)

## INGROWN NAIL SPLITTING SCISSORS

Product Number: 97-0555

6" Ingrown Nail Splitting Scissors One Pointed Blade One Serrated Blade

Sklar's® Ingrown Nail Splitting Scissors are a podiatric instrument used in procedures for the treatment of ingrown nails. These non-ratcheted, finger ring scissors have one pointed, sharp tip and one rounded, blunt tip. The blade of the pointed tip is serrated, while the blade of the blunt tip is smooth. This product is primarily used for cutting wedges in the nail plate as part of removing the ingrown section. They are straight and the overall length is 6 inches.



### Additional Information

Length	6"
Tip Configuration	Sharp/Blunt
Instrument Type	Scissors
Curvature	Straight
Working Surface Style	Serrated
Material	Stainless Steel
Disposable or Reusable	Reusable
Sterile or Non-Sterile	Non-Sterile
Latex or Latex-Free	Latex-Free
Grade	Premium OR-Grade

© 2016 Sklar Surgical Instruments. All Rights Reserved.

889 South Matlack Street, West Chester, PA 19382 • Phone: 610.756.7865 • Fax: 610.430.3941 • [Service@SklarCorp.com](mailto:Service@SklarCorp.com)